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Care-Lock® priključak kompatibilan je samo s Care-Lock® kompletom za prijenos.

Ulaz za šiljak kompatibilan je sa svim Fresenius Kabi kompletima za prijenos. Mogu se koristiti i drugi sustavi za probadanje, a korisnik mora testirati njihovu prikladnost. Na temelju testova provedenih s drugim sustavima za prijenos na tržištu, poznato je da sljedeći kompleti nisu kompatibilni s ulazom za šiljak vrećice Careflex® Duo: Pump set for ARTHREX-System tvrtke Meise Medizintechnik (Ref.: AR 6410), Pump Tube for fms solo-duo tvrtke Meise Medizintechnik (Ref.: AR 7410), Tubing for Cystoscope or Resectoscope tvrtke Polfa-Lublin, Normoflo Irrigation Warming Set tvrtke Smiths Medical (Ref.: IRI-600), AHTO Disposable Suction Irrigator Tube Set with Tip tvrtke Stryker (Ref.: 250-070-620), Flocontrol Tube Set Pump SIL, Reusable tvrtke Stryker (Ref.: Z1428-39), TUR Irrigation tvrtke Techno Med (Ref.: TY30002SR), Fluidio Irrigation URO Set tvrtke The Surgical Company International B.V. (Ref.: 670800-B), Fluidio Irrigation URO Spike Set tvrtke The Surgical Company International B.V. (Ref 680803-B), Single Use Tubing Set for Pump PV-5201 tvrtke Vimex Sp.Z.o.o. (Ref.: PV-5201SUTS).

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INSTRUCTION FOR USE - READ CAREFULLY!

Fresenius Kabi AG
61346 Bad Homburg
Germany

NaCl 0.9% Fresenius

Solution for irrigation

Composition	1 Litre contains:	
	Sodium chloride	9 g
	Hydrochloric acid (pH adjustment)	
	Sodium hydroxide (pH adjustment)	
	Water for injections in bulk, Ph. Eur.	
	Equivalent to:	Na ⁺ 154 mmol/l Cl ⁻ 154 mmol/l
	pH value:	5.0 - 7.0
	Titration acidity	< 1 mmol NaOH/l
	Theoretical osmolarity	308 mosm/l

Indications	Plasmaisotonic (physiological) saline solution for internal and external application such as
	<ul style="list-style-type: none">- Postoperative bladder irrigation after all urological interventions- Irrigation of the gastrointestinal tract- Irrigation of fistulas and drainages- Wound treatment and to moisten cloths and dressings.- Perfusion of the extracorporeal system in hemodialysis

Contraindications	Hypernatremia, hyperhydration. Caution in case of hyperchloreaemia.
--------------------------	--

Precautions	Use with caution in patients with hyperhydration, hypernatremia, and hyperchloreaemia.
--------------------	--

Side-effects

Hyperhydration, hypernatremia, hyperchloremia and metabolic acidosis may occur in case of absorption of large amounts of NaCl 0.9% Fresenius.

Interaction with other drugs

None known.

Dosage

Dosage and quantity to be applied will result from the type of intervention.

Description and connectivity

The Careflex® Duo bag is a polyolefin bag with a spiking port and the Care-Lock® port.

The Care-Lock® port is only compatible with the Care-Lock® transfer set.

The spiking port is compatible with all Fresenius Kabi transfer sets. Other spiking systems may be used and have to be tested for suitability by the user. Based on tests conducted with other transfer systems on the market the following sets are known to be not compatible with the spiking port of the Careflex® Duo bag: Pump set for ARTHREX-System by Meise Medizintechnik (Ref: AR 6410), Pump Tube for fms solo-duo by Meise Medizintechnik (Ref: AR 7410), Tubing for Cystoscope or Resectoscope by Polfa-Lublin, Normoflo Irrigation Warming Set by Smiths Medical (Ref: IRI-600), AHTO Disposable Suction Irrigator Tube Set with Tip by Stryker (Ref: 250-070-620), Flocontrol Tube Set Pump SIL, Reusable by Stryker (Ref: Z1428-39), TUR Irrigation by Techno Med (Ref: TY30002SR), Fluido Irrigation URO Set by The Surgical Company International B.V. (Ref: 670800-B), Fluido Irrigation URO Spike Set by The Surgical Company International B.V. (Ref 680803-B), Single Use Tubing Set for Pump PV-5201 by Vimex Sp.Z.o.o. (Ref: PV-5201SUTS).

Use of non-compatible transfer systems may lead to leakages or loss of spike during use and impossibility of irrigation.

An instruction how to spike the bag can be found at the end of this instruction for use.

Presentation

Polyolefin bags (Careflex® Duo) 5 x 2,000 ml
4 x 3,000 ml
2 x 5,000 ml

Additional Information

Not for intravenous use or injection!

The sterility of NaCl 0.9% Fresenius is guaranteed only until the package is opened. The solution must be used immediately after opening.

Use only if the solution is clear and free from particles and if the container is undamaged.

For single use only.

Any unused portion of the solution must be discarded.

Re-use of medical single use equipment involves a risk of contamination of the equipment. This may cause infection in the patient or the user, which may lead to illness or death.

Keep out of the sight and reach of children!

Do not use NaCl 0.9% Fresenius after the expiry date.

Date of revision

September 2021

	CE marking of conformity with Medical Devices Directive 93/42/EEC (as amended) with the intervention of Notified Body 0123
	Indicates the presence of a sterile fluid path within the medical device in cases when other parts of the medical device, including the exterior, might not be supplied sterile.
	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.

	Indicates the medical device manufacturer.
	Indicates the need for the user to consult the instructions for use.
	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Indicates the date after which the medical device is not to be used.
	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Indicates a medical device that should not be used if the package has been damaged or opened.
	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	This product is not for intravenous use or injection.
	Indicates a medical device that needs protection from light sources.
	Indicates a medical device that needs to be protected from moisture.
	Indicates the temperature limits to which the medical device can be safely exposed. Store at room temperature.


Spiking the bag

For spiking, the bag can be hung on an IV pole or alternatively placed on a flat surface. Administer the content in the hanging position.




Check product name, bag integrity and the expiry date. Inspect the container for damage or leakage and for particles. **Do not use if the container is damaged or if the solution contains particles.**

Identify the blue spiking port.

Break off the blue tamper-evident cover from the port.
Do not touch the opened surface!

Close clamp of transfer set if applicable. Insert the spike **centered** and **straight** until the blue plastic collar of the port meets the shoulder of the spike.

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