# GLUCOFAST® DUO

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PZN 15405027

PZN 15401928

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
## BLOOD GLUCOSE METER

## OWNER'S MANUAL

For self-testing and professional use.

## Dear Glucofast® Duo System Owner:

Thank you for purchasing the Glucofast® Duo Blood Glucose Monitoring System. This manual provides important information to help you to use the system properly. Before using this product, please read the following contents thoroughly and carefully.

Regular monitoring of your blood glucose levels can help you and your doctor gain better control of your diabetes. Due to its compact size and easy operation, you can use the Glucofast® Duo Blood Glucose Monitoring System to easily monitor your blood glucose levels by yourself anywhere, any time.

If you have other questions regarding this product, please contact the local customer service or place of purchase.

## Intended Use

This system is intended for use outside the body (in vitro diagnostic use) by people with diabetes at home and by health care professionals in clinical settings as an aid to monitoring the effectiveness of diabetes control. It is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples from the fingertips, and from venous, arterial and neonatal blood sample. It should not be used for the diagnosis of or screening for diabetes mellitus.

The test for hematocrit (HCT) as part of the system, is intended for use in the in vitro quantification of packed red blood cell volume fraction in capillary whole blood as an aid in monitoring the status of total volume of red blood cells. The test reading of hematocrit (HCT) is used only to determine whether the blood test sample is within the acceptable range of Blood Glucose Monitoring System. It should not be used for the diagnosis of anaemia or erythrocytosis.

## Test Principle

Your system measures the amount of glucose in whole blood. The glucose testing is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current, calculates the glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

## **IMPORTANT SAFETY INSTRUCTIONS**

### **Read before the first use.**

1. Use this device **ONLY** for the intended use described in this manual.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or if it is damaged.
4. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only. Always consult your doctor to have the results interpreted.
5. The blood glucose test strip can be used for the testing of newborns.
6. Before using this device to test blood glucose, read all instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
7. Keep the device and testing equipment away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
8. Use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.
9. Do **NOT** use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with accurate operation.
10. Proper maintenance and periodic control solution testing are essential to the longevity of your device. If you are concerned about your accuracy of measurement, please contact the local customer service or place of purchase for help.

**Keep these instructions in a safe place!**

## TABLE OF CONTENTS

<b>BEFORE YOU BEGIN.....</b>	<b>38</b>
Important Information	38
Meter Overview	39
Test Strip	41
<hr/>	
<b>SETTING THE METER.....</b>	<b>42</b>
<hr/>	
<b>THE THREE MEASURING MODES.....</b>	<b>44</b>
<hr/>	
<b>QUALITY CONTROL TESTING.....</b>	<b>45</b>
When should the Control Solution Test be Performed?	45
Performing a Control Solution Test	45
<hr/>	
<b>TESTING WITH BLOOD SAMPLE.....</b>	<b>47</b>
Preparing the Lancing Device for Blood Testing	47
Preparing the Puncture Site	47
Performing a Blood Glucose Test	48
<hr/>	
<b>METER MEMORY .....</b>	<b>50</b>
Reviewing Test Results	50
Reviewing Blood Glucose Day Average Results	50
<hr/>	
<b>DOWNLOADING RESULTS ONTO A COMPUTER .....</b>	<b>51</b>
Data Transmission via Cable	51

**MAINTENANCE.....52**

Battery	52
Replacing the Battery	52
Caring for your Meter	53
Caring for your Test Strips	54
Important Control Solution Information	55

---

**MEASUREMENT RESULT READINGS.....55**

Result Readings	55
Reference Values	56

---

**SYSTEM TROUBLE SHOOTING.....56**

Error Messages	56
Troubleshooting	57

---

**SYMBOL INFORMATION.....59**

---

**SPECIFICATIONS.....60**

## **BEFORE YOU BEGIN**

### **Important Information**

Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.

If your blood glucose results are lower or higher than usual, and you do not have any symptoms of illness, first repeat the test. If you have symptoms or continue to get results which are higher or lower than usual, follow the treatment advice of your healthcare professional.

Use only fresh whole blood samples to test your blood glucose. Using other substances will lead to inaccurate results.

If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all the instructions given in this owner's manual, contact your healthcare professional.

We do not recommend using this product on severely hypotensive individuals or patients in shock. Please consult your healthcare professional before use.

The measurement unit used for indicating the concentration of blood or plasma glucose will have molarity (mmol/L). The approximate calculation rule for conversion of mmol/L in mg/dL is:

mg/dL	Divided by 18	= mmol/L
mmol/L	Times 18	= mg/dL


#### **For example:**

- 1)  $120\text{mg/dL} \div 18 = 6.6\text{mmol/L}$
- 2)  $7.2\text{mmol/L} \times 18 = 129\text{mg/dL}$  approximately.

The HCT is the volume percentage (%) of red blood cells in blood sample.

(Note: Hematocrit = HCT)

## Meter Overview

### 1. Test Strip Slot

Insert test strip here to turn the meter on for testing.

### 2. Display Screen

### 3. Main Button (M)

### 4. Test Strip Ejector

Eject the used strip by pushing up this button.

### 5. Data Port


Download test results with a cable connection.

### 6. SET Button (S)

Enter and confirm the meter settings.

### 7. Battery Compartment

## Display Screen

1. Glucose Symbol

2. Test Result

3. Error Message

4. Low Battery Symbol

5. Day Average

6. Memory Mode

7. Date/Time

8. HCT Level

9. Measurement Unit

10. Alarm Reminder

11. Auto QC Mode/QC Mode  
QC- control solution test

**Measuring Mode**

AC – before meal

PC – after meal

Gen – any time of day


12. Buzzer Symbol

13. Test Strip Symbol

14. Blood Drop Symbol

## Test Strip

### 1. Absorbent Hole

Hold the drop of blood on this field. The blood is automatically moved in slowly.

### 2. Confirmation Window

Provides if sufficient blood has been sucked.

### 3. Test Strip Handle


At this point, the test strip can be touched while inserting it.

### 4. Contact bars


Insert this end of the test strip into the meter and push it in until it stops.

## ATTENTION

**The front side of test strip should face up when inserting test strip.**

Test results might be wrong if the contact bar is not fully inserted into the test slot.

## !! NOTE


The Glucofast® Duo monitor should only be used with Glucofast® Duo Test Strips. Using other test strips with this meter can produce inaccurate results.

## SETTING THE METER

Before you use your device for the first time, or if you change the battery, you should check these settings and update them if necessary. Make sure you save the settings you want after following the steps below.

### Entering the Setting Mode

Start with the meter off (no test strip inserted). **Press S.**

#### 1. Setting the date

The sequence of the date setting is: YEAR -> MONTH -> DAY. With the YEAR / MONTH / DAY flashing in sequence, **press M** to select the correct number. **Press S.**

#### 2. Setting the time format

**Press M** to select the desired time format - 12h or 24h. **Press S.**

#### 3. Setting the time

With the HOUR/MINUTE flashing in sequence, **press M** until the correct hour/minute appears. **Press S.**

#### 4. Setting the buzzer

With the buzzer displays, **press M** to switch between "On" and "OFF". **Press S.**

#### 5. Deleting the memory

With "dEL" and a **M** on the display, **press M** and select "no" to keep the results in memory then **press S** to skip. To delete all the results, **press M** and select "YES" to delete all the memory records.

#### 6. Setting the reminder alarm

Your meter has four reminder alarms. The meter will display "OFF" and "AL1". If you don't want to set an alarm, **press S** to skip this step; or **press M** to select "On", then **press S.**

With the hour/minute flashing in sequence, **press M** to select the correct hour/minute. **Press S** and go to the next alarm setting.

## NOTE:

When the alarm beeps, press M to switch it off. Otherwise, it will beep for 2 minutes then switch off.

## 7. Setting the backlight

The default setting for meter backlight ("BL") is set to "On". **Press M** to switch between "On" and "OFF". **Press S**.

**Congratulations! You have completed all settings!**

## NOTE:

- These parameters can **ONLY** be changed in the setting mode.
- If the meter is idle for 3 minutes during the setting mode, it will switch off automatically.

## THE THREE MEASURING MODES

The meter provides you with three modes for measuring.

MODE	TIME TO USE
General (Displayed as "Gen")	At any time of the day, without consideration on the last meal
AC 	Before eating, no meal since at least 8 hours.
PC 	2 hours after a meal.

**You can choose between the different modes:**

1. Start with the meter switched off. Insert a test strip to turn on the meter. The screen will display a flashing "🔥" and "Gen".
2. **Press M** to switch between General "Gen", "AC ", "PC " mode.

## QUALITY CONTROL TESTING

### When Should the Control Solution Test be Performed?

Our Control Solution contains a known amount of glucose that reacts with test strips and is used to ensure your meter and test strips are working together correctly.


#### When Should the control solution test be performed?

- if it is mandatory following the local regulations in your country,
- if you suspect your meter or test strips are not working properly,
- if your blood glucose test results are not consistent with how you feel, or if you think the results are not accurate,
- to practice the testing process, or
- if you have dropped or think you may have damaged your meter.

Test strips, control solution, lancing device or sterile lancets may not be included in the kit (please check the contents in your product box). They can be purchased separately. Please make sure you have all items needed for a blood glucose test beforehand.

### Performing a Control Solution Test

To perform a control solution test, you will need: the **meter**, a **test strip** and a **control solution vial**.

#### 1. Insert the test strip to turn on the meter


Wait for the meter to display "C" and "▲".

#### 2. Apply control solution

Shake the control solution vial thoroughly before use (**picture 1**).

Squeeze out the first drop and wipe it off, then squeeze out another drop and place it on the tip of the vial cap (**picture 2**).

Hold the meter to move the absorbent hole of the test strip to touch the drop (**picture 3**). Once the confirmation window fills completely, the meter will begin counting down.

**!! NOTE:**

- To avoid contaminating the control solution, do not directly apply control solution onto a strip.
- Your device will tag this measurement as "QC" test automatically.

**3. Read and compare the result**

After counting down to 0, the control solution test result will appear on the display. Compare this result with the range printed on the test strip vial and it should fall within this range. If not, please read the instructions again and repeat the control solution test.

With "QC" displayed, the meter will store your test result in memory under "QC".

**!! NOTE:**

- The control solution range printed on the test strip vial is for control solution use only. It is not a recommended range for your blood test level.
- See the **MAINTENANCE section** for important information about your control solutions.

## TESTING WITH BLOOD SAMPLE

### !!WARNING:

To reduce the chance of infection:

- Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.

### Preparing the Lancing Device for Blood Testing

Please follow the instructions in the lancing device insert for collecting a blood sample.

### Preparing the Puncture Site

Stimulating blood perfusion by rubbing the puncture site before blood extraction has a significant influence on the accuracy of the glucose value obtained.

Blood from a site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. If the puncture site was rubbed prior to blood extraction, the difference is significantly reduced.

#### **Please follow the suggestions below before obtaining a drop of blood:**

- Wash and dry your hands before starting.
- Select the puncture site at fingertips.
- Rub the puncture site for about 20 seconds before puncture.
- Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry**.
- **Fingertip testing**

Press the lancing device's tip firmly against the lower side of your fingertip. Press the release button to prick your finger, then a click indicates that the puncture is complete.


## !! NOTE:

Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.

### Performing a Blood Glucose Test

To perform a blood glucose test, you will need: the **meter**, one blood glucose **test strip**, the **lancing device** and **lancets**.


#### 1. Insert the test strip to turn on the meter

Wait for the meter to display the "E" and "H", and "GLU".

#### 2. Select the appropriate measuring mode by pressing M.

#### 3. Obtaining a blood sample

Use the pre-set lancing device to puncture the desired site. The size of the drop should be at least as big as ● (actual size), which is 1.0 microliter ( $\mu\text{L}$ ) of volume. Gently squeeze the punctured area to obtain another drop of blood. Be careful **NOT** to smear the blood sample.

#### 4. Apply the sample

Gently apply the drop of blood to the absorbent hole of the test strip at a tilted angle. Confirmation window should be completely filled if enough blood samples have been applied. Do **NOT** remove your finger until you hear a "beep" sound.

**NOTE:**


- Do not press the punctured site against the test strip or try to smear the blood.
- If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to start a new test.
- The confirmation window should be filled with blood before the meter begins to count down. **NEVER try to add more blood to the test strip after the drop of blood has moved away. Discard the used test strip and retest with a new one.**
- If you have trouble filling the confirmation window, please contact your health care professional or the local customer service for assistance.

### 5. Read Your result

The result of the blood glucose test will appear after the meter counts down to 0. The result will be stored in your memory automatically together with the HCT value.

### 6. Remove the used test strip and lancet

Eject your test strip by pushing the eject button on the side. Use a sharp bin to dispose of used test strips. The meter will switch itself off automatically.

**Always follow the instructions in the lancing device insert when removing the lancet.**

**WARNING:**

The used lancet and test strip may be biohazardous. Please discard them carefully according to your local regulations.

## METER MEMORY

The Glucofast® Duo stores the 1000 most recent blood glucose and HCT test results, dates and times in its memory. To enter the meter memory, **start with the meter switched off.**

### Reviewing Test Results

#### 1. Press and release M.

"**M**" will appear on the display. **Press M** again, and the first reading you see is the last blood glucose and HCT result along with date, time and the measuring mode.

#### 2. Press M to recall the test results stored in the meter each time you press. **After the last glucose test result, press M again and the meter will switch off.**

### Reviewing Blood Glucose Day Average Results

1. **Press and release M.** When "**M**" appears on the display, **keep pressing M** for 3 seconds until the flashing "**DAY AVG**" appears. **Release M** and then your 7-day average result measured in general mode will appear on the display.

2. **Press M to review** 14-, 21-, 28-, 60- and 90- day average results stored in each measuring mode in the order of Gen, AC, and then PC.

3. **Exit the meter memory. Keep pressing M** and the meter will switch itself off after displaying the last test result.

#### !! NOTE:

- Any time you wish to exit the memory, **keep pressing M** for 3 seconds or leave it without any action for 3 minutes. The meter will switch off automatically.
- Control solution results are **NOT** included in the day average.

# **DOWNLOADING RESULTS ONTO A COMPUTER**

## **Data Transmission via Cable**

You can use the meter with a USB cable (optional) and the Healthcare Software System to view test results on your personal computer. To learn more about the Healthcare Software System or to obtain a USB cable separately, please contact local customer services or the place of purchase for assistance.

### **1. Obtaining the required cable and installing the software**

You can purchase a standard USB cable commercially. The free Health Care Software can be found on the TaiDoc website ([www.taidoc.com](http://www.taidoc.com)).

### **2. Connecting to a personal computer**

Connect the cable to a USB port on your computer with the meter switched off; connect the other end of the USB cable to the meter data port. "PC" will appear on the meter display, indicating that the meter is in communication mode.

### **3. Data transmission**

To transmit data, follow the instructions provided with the software. Results will be transmitted with date and time. Remove the cable and the meter will automatically switch off.

#### **!!! WARNING:**

While the meter is connecting to the PC, it will be unable to perform a test.

## MAINTENANCE

### Battery

Your meter comes with two 1.5V AAA size alkaline batteries.

#### Low Battery Signal

The meter will display one of the messages below to alert you when the meter power is getting low.

1. The "🔋" symbol appears. The meter is functional and the result remains accurate, but it is time to change the batteries.
2. The "🔋" symbol appears with "E-b" and "⚠️": The power is not enough to do a test. Please change the batteries immediately.

#### Replacing the Battery

**To replace the batteries, make sure the meter is turned off.**

1. Press the edge of the battery cover and lift it up to remove.
2. Remove the old batteries and replace with two 1.5V AAA size alkaline batteries.
3. Close the battery cover. If the batteries are inserted correctly, you will hear a "beep" afterwards.

## !!! NOTE:

- Replacing the batteries does not affect the test results stored in the memory.
- As with all small batteries, these batteries should be kept away from children. If swallowed, promptly seek medical assistance.
- Batteries might leak chemicals if unused for a long time. Remove the batteries if you are not going to use the device for an extended period (i.e., 3 months or more).
- Properly dispose of the batteries according to your local environmental regulations.

## Caring for Your Meter

To avoid contaminating the meter and test strips, wash your hands before each use and dry them thoroughly.

### Cleaning

1. To clean the meter externally, wipe it with a slightly damp cloth. Use tap water or a mild detergent. Then dry the device with a soft and dry cloth. **DO NOT** hold the device under running water.
2. Do **NOT** use organic solvents to clean the meter.

### Meter Storage

- Storage conditions: -20°C to 60°C, between 10% and 93% relative humidity (non-condensing).
- Always store or transport the meter in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

### Meter Disposal

A used meter should be treated as contaminated and may carry a risk of infection during measurement. The batteries in this used meter should be removed and the meter should be disposed in accordance with local regulations.

The meter falls inside the scope of the European Directive 2002/96/ EC-Directive on waste electrical and electronic equipment (WEEE). Do not discard in household waste.

## Caring for Your Test Strips

- Storage conditions: 2°C to 30°C, between 10% and 85% relative humidity (non-condensing). **Do NOT freeze.**
- Store your test strips in their original vial only. Do not transfer to another container.
- Store test strip packages in a cool dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands. Use each test strip immediately after removing it from the vial.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on the vial label when you first opened it. Discard remaining blood glucose test strips after 6 months.
- Do not use test strips after the expiration date. This may cause inaccurate results.
- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.


**For further information, please refer to the test strip package insert.**

## Important Control Solution Information

- Use only Glucofast<sup>+</sup> Duo control solutions with your meter.
- Do not use the control solution beyond the expiration date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months
- It is recommended that the control solution test be done at room temperature 20°C to 25°C. Make sure your control solution, meter and test strips are at this specified temperature range before testing.
- Shake the vial before use, discard the first drop of control solution and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 2°C to 30°C. Do NOT freeze.

## MEASUREMENT RESULT READINGS

### Result Readings

MESSAGE	WHAT IT MEANS
	< 10mg/dL (< 0.56mmol/L)
	between 10mg/dL - 239mg/dL (0.56mmol/L - 13.2mmol/L)
	≥ 240mg/dL (≥ 13.3mmol/L)
	> 600mg/dL (> 33.3mmol/L)

## Reference values

The meter provides you with plasma equivalent results.

TIME OF DAY	NORMAL PLASMA GLUCOSE RANGE FOR PEOPLE WITHOUT DIABETES
Fasting and before meal	<100mg/dL (5.6mmol/L)
2 hours after meals	<140mg/dL (7.8mmol/L)

Source: American Diabetes Association. Standards of Medical Care in Diabetes- 2018 Jan; 41 (Supplement 1): S1-S2.

**Please consult your doctor to determine a target range that works best for you.**

## SYSTEM TROUBLESHOOTING

### Error Messages

If you follow the recommended action but the problem persists, please contact your local customer service.

**Never try to disassemble the meter or repair it yourself.**

MESSAGE	WHAT IT MEANS	WHAT TO DO
<b>E-b</b>	Appears when the battery is too low	Replace the battery immediately
<b>E-U</b>	Appears when a used test strip is inserted	Repeat with a new test strip
<b>E-t</b>	Appears when ambient temperature is above or below system operation range	System operation range is 10°C to 40°C. Repeat the test after the meter and test strip are in the above temperature range

<b>E-F</b>	Appears when test strip is removed while counting down, or insufficient blood volume	Review the instructions and repeat test with a new strip. If the problem persists, please contact the local customer service for help
<b>E-O</b> <b>E-A</b> <b>E-E</b> <b>E-C</b>	Problem with the meter	Repeat the test with a new test strip. If the meter still does not work, please contact the customer service for assistance

## Troubleshooting

### 1. If the meter does not display a message after inserting a test strip:

POSSIBLE CAUSE	WHAT TO DO
Batteries exhausted	Replace the batteries
Test strip inserted upside down or incompletely	Insert the test strip with contact bars end first and facing up
Defective meter or test strips	Please contact customer services

## 2. If the test does not start after applying the sample:

POSSIBLE CAUSE	WHAT TO DO
Insufficient blood sample	Repeat the test using a new test strip with larger volume of blood sample
Defective test strip	Repeat the test with a new test strip
Sample applied after automatic switch-off (3 minutes after last user action)	Apply sample only when flashing "⚡" appears on the display
Defective meter	Please contact customer services

## 3. If the control solution testing result is out of range:

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test	Read instructions thoroughly and repeat the test again
Control solution vial was poorly shaken	Shake the control solution vigorously and repeat the test again
Expired or contaminated control solution	Check the expiry date of the control solution
Control solution that is too warm or too cold	Control solution, meter, and test strips should be at room temperature 20°C to 25°C before testing
Defective test strip	Repeat the test with a new test strip
Meter malfunction	Please contact customer services
Improper working of meter and test strip	Please contact customer services

## SYMBOL INFORMATION

SYMBOL	REFERENT
	<i>In vitro</i> diagnostic medical device
	Consult instructions for use
	Temperature limit
	Use-by date
	Batch code
	Serial number
	Manufacturer
	Authorized representative in the European Community
	Caution
	Humidity limitation
	CE mark
	Catalogue number
	Do not discard in household waste

## SPECIFICATIONS

<b>Model No.:</b>	TD-4279
<b>Dimension &amp; Weight:</b>	96 (L) x 61 (W) x 26 (H) mm, 67.2g
<b>Power Source:</b>	Two 1.5V AAA alkaline batteries
<b>Display:</b>	LCD with backlight
<b>Memory:</b>	1000 measurement results with respective date and time
<b>External Output:</b>	USB

Auto electrode insertion detection

Auto sample loading detection

Auto reaction time count-down

Auto detection control solution mode

Auto switch-off after 3 minutes without action

Temperature Warning

One code system

**Operating Condition:** 10°C to 40°C,  
between 10% and 85% R.H. (non-condensing)

**Meter Storage/Transportation Conditions:** -20°C to 60°C,  
between 10% and 93% R.H. (non-condensing)

**Strip Storage/Transportation Conditions:** 2°C to 30°C  
and 10% to 85% R.H. (non-condensing)

**Measurement Units:** default mg/dL or mmol/L

**Measurement Range:** 10 to 600mg/dL (0.56 to 33.3mmol/L)

**Expected Service Life:** 5 years

**Operating Altitude:** Up to 2000m, for indoor use

**Degree of Pollution:** Pollution degree 2

This device has been tested to meet the electrical and safety requirements of:  
IEC/EN 61010-1, IEC/EN 61010-2-101, IEC/EN 61326-1, IEC/EN 61326-2-6.

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Zur Eigenanwendung.  
For self-testing.

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**CARDIMAC GmbH**
**MedNet GmbH**
**TaiDoc Technology Corporation**