- Ⓓ Gebrauchsanweisung
- Ⓔ Instructions for use
- Ⓕ Mode d'emploi
- Ⓖ Instrucciones para el uso
- Ⓘ Istruzioni per l'uso
- Ⓝ Gebruiksaanwijzing

- Ⓟ Instruções de utilização
- Ⓡ Brugsanvisning
- Ⓢ Bruksanvisning
- Ⓛ Instrukcja obsługi
- Ⓜ Οδηγίες χρήσης
- Ⓒ Návod k použití

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## ☉ Hot and cold pillow with loops

The pillow, with a filling of natural lavender clay beads, can help relax and reduce stress and can be used anywhere on your body thanks to its practical loops. The pillow is suitable for hot and cold application. Please make sure to observe the following information.

**Heat application in the microwave:** Since the actual performance of various microwave oven manufacturers with the same specification is different, we recommend testing the optimal warm-up time.

Place the pillow flat into the microwave and distribute the filling evenly. Heat the pillow at max. 600 watts for about 2 minutes. If the pillow is not warm enough, knead it and heat it for another 30 seconds. Repeat the procedure in steps of 30 seconds until the desired temperature is reached. Kneading is important to distribute the heat evenly and to ensure that there is no overheating of individual spots.

### Cold application in the refrigerator:

Place the neck and shoulder pillow in the fridge or freezer for at least 2 hours before use. We recommend to put the pillow into a bag before cooling to protect it against moisture and odours.

**Cleaning:** If the pillow requires cleaning, it can be wiped with a damp cloth.

### IMPORTANT INFORMATION:

**Caution:** Do not heat in the oven or on the grill. Never leave the pillow unattended during heating.

Do not exceed the specified wattage.

**Risk of burns:** Check the temperature of the heat pillow BEFORE you use it on yourself or others. The pillow is not suitable for children.

**Risk of fire: in case of overheating and drying out.** The turntable of the microwave must be clean and move freely.

Follow the instructions for use exactly. Keep away from children.

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