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It is generally assumed that different types of wounds are colonized by microorganisms. This does not cause a disease, but decon- tamination is always desirable in order to prevent wound infections and to provide optimal support for the wound healing process by cleaning the wound thoroughly and by sparing of viable tissue.

General information – ActiMaris®Forte Wound rinsing solution (300ml)

for the decontamination, cleaning and moistening of acute and chronic wounds, skin, mucous membranes and first and second degree burns.

1. Product profile

ActiMaris® Forte – Wound rinsing solution consists of sea salt, water and reactive oxygen.

Reactive oxygen eliminates wound typical pathogens, without damaging the tissue. The oxidative effect of singlet oxygen (102) is also responsible for the neutralization of wound odours.

The alkaline pH - value of ActiMaris® Forte – Wound rinsing solution additionally supports the controlled release of reactive oxygen, by providing the OH groups (as part of the naturally balanced redox system).

Sea salt contains sodium chloride as its main ingredient and over 80 of the other different components of sea water. In addition to potassium, calcium, and magnesium, sulphur, bromine, iron and essential traces of iodine, zinc, manganese and selenium are also present. As in seawater, these components are present in their natural complex compounds and in physiologically correct concentrations so they favourably influence the cell milieu and (wound) metabolism.

The increased salinity of ActiMaris® - Wound rinsing solution allows for the rinsing of wound oedema so that swelling subsides rapidly.

2. Fields of application

ActiMaris® Forte - Wound rinsing solution is suitable:

For cleaning, rinsing, moistening and decontamination of

- acute mechanical wounds (skin lacerations, bites, cuts, abrasions, tears and bruises) and postoperative wounds.
- chronic wounds (decubitus ulcers, arterial-venous ulcers, diabetic ulcers)
- necrotic, malodorous wounds and tumor ulcer
- thermal and chemical wounds (first and second degree burns)
- entry portals of urological catheters and PEG tubes and drains.

Intraoperative cleaning and rinsing of wounds.

Moistening of bandages and wound dressings.

During dressing changes to loosen encrusted bandages or other encrusted wound dressings.

ActiMaris® Forte – Wound rinsing solution is suitable particularly for the initial cleaning of severe fibrinous or purulent wounds, as well as heavily infected or strong foul-smelling wounds.

3. Instructions for use

ActiMaris® Forte – Wound Rinsing Solution can be applied at room temperature. During washing of wounds with increased pain inclination, it should be applied at body temperature (rinsing solution is to be heated to approximately 37°C).

Before using ActiMaris® Forte – Wound rinsing solution, the wounds should first be rinsed and cleaned with ActiMaris® Sensitiv – Wound rinsing solution, in order to remove dirt as well as easily removable coatings.

The subsequent deep cleansing and odour reduction will be with ActiMaris® Forte –Wound rinsing solution; for this purpose the soaked sterile bandages or wound dressings are applied. (see also section 8: Interactions).

It should be applied primarily as classic wet/dry phase (according to G. Kammerlander).

The recommended application time of the wet phase in case of severe fibrinous or purulent, infected or strong foul-smelling wounds is → 5-10 minutes.

The duration and frequency of application should be adapted to the individual wound conditions. The use of the wound rinsing so- lution is not limited in time and frequency. Necrosis and coatings can be removed mechanically, depending on the wound condition. Primary surgical sanitation is advised with thick, adherent, fibrinous or necrotic wound coatings. After the improvement of the wound condition, it is advised to continue further wound cleansing with ActiMaris® Sensitiv – Wound rinsing solution until the wound is optically clean.

4. Tissue tolerance and biocompatibility

ActiMaris® Forte – Wound rinsing solution is suitable for use even on sensitive / irritated skin and mucosa due to its non-irritating and well-tolerated application. There is no Inhibition of granulation and epithelialisation. Dermatologically approved.

5. Side effects

No side effects known at present.

6. Contraindications

ActiMaris® Forte – wound rinsing solution is not to be applied in case of known or suspected allergy to any of the ingredients.

7. Restrictions on use:

Pregnancy and lactation:

There are no indications of mutagenic or embryo toxic effects of the ingredients. Due to the lack of studies and long-term clinical experience with pregnant and breast feeding women, ActiMaris® Forte – Wound rinsing solution should be used in these cases only after an accurate medical assessment.

Infants and newborns:

Due to the lack of studies and long-term clinical experience with pregnant and breast feeding women,, ActiMaris® Forte – Wound rinsing solution should be used in these cases only after an accurate medical assessment.

8. Interactions

ActiMaris® Forte – Wound rinsing solution works due to the oxidation of the reactive oxygen's content.

To avoid possible interactions, the rinsing solution should not be mixed with or used in combination with substances containing silver, biguanides, taurolidine or iodine-containing solutions and wound dressings. This shall also apply even to wetting of bandages and wound dressings; here only sterile, not active bandaging materials such as e.g. active agent free dressings, gauze, pads, sponges, gels, hydro fibres, alginates, hydrocolloids and similar should be used. Simultaneous application of the rinsing solution with other ActiMaris® - products (Wound gel and Sensitiv – Wound rinsing solution) is possible.

If ActiMaris® Forte – Wound rinsing solution is used as the classic wet/dry phase and the corresponding drying time is taken into account, there are no interactions with consecutively applied wound care products

9. General safety instructions

ActiMaris® Forte – Wound rinsing solution is suitable only for external wounds on skin or mucosa. Do not use for infusion or injection! Do not swallow!

Use only closed and undamaged containers.

Keep out of reach of children!

It is possible for coloured fabrics to fade after contact with ActiMaris® - products.

10. Composition and technical information:

Sea salt, water.

Clear, colorless, basic and hypertonic sea salt solution with characteristic smell; contains 3% NaCl (Sal Maris) and 0,2% NaOCl.

11. Storage specification, duration, expiry date

ActiMaris® Forte – Wound rinsing solution has to be stored at room temperature between +5° and +25°C. Protect the bottle from the direct sunlight.

Shelf life is 30 months; Expiry date on the packaging.

Period of use is 3 months after opening the bottle. Close carefully after use.

ActiMaris AG

simple medical solutions

Manufacturer

ActiMaris AG
Technikumstrasse 14
CH-9470 Buchs
admin@actimaris.com
Tel.: +41 71 505 75 25

Lot number
Order number
Follow instructions for use
Use by: year and month
Storage (temperature)
Protect from direct sunlight
Manufacturer

Date of information: 06/2013 - V3.0
Medical product: Klasse II b, CE 0408