Starlyf® ULTRASONIC

ENGLISH

Thank you for purchasing the Starlyf® Ultrasonic, the tooth cleaner that uses vibration to gently remove stains, plaque and tartar. Before using the product for the first time, please read all these instructions to ensure you use the product correctly.

Contents:
1. Warnings
2. The Set Contains
3. Instructions for Use
4. Cleaning and Maintenance
5. Quality Warranty

1. WARNINGS

Always read the instructions before using the Starlyf® Ultrasonic.

. It runs on 1 AAA battery (not included).

. If the battery is in bad condition, do not use it.

Do not use the product if it is damaged or defective.

Do not use the Starlyf® Ultrasonic if your hands are wet.

Always keep unused batteries in the original packaging and away from metal objects that could cause them to short-circuit.

Remove the battery from the device if you do not plan to use it for an extended period of time.

Never dispose of batteries in fire, as they may explode.

Do not recharge the battery.

Do not tamper with the device in any way. The device should be handled or repaired only by an authorized professional.

. This product is not a toy and should be kept out of reach of children.

This product can be used by children aged eight and over, people with reduced physical, sensory or mental abilities, and people without experience or knowledge, provided they are supervised or given instructions on how to use the device safely and understand the hazards involved in its use. Children should not be allowed to play with the device. Children should not clean or do maintenance on the device without adult supervision. Consult your dentist before using the Starlyf® Ultrasonic on children under eight years of age.

Do not use the product for any purpose not recommended in this manual.

2. The Set Contains

- Starlyf® Ultrasonic unit
- Stainless steel head
- Silicone head
- Tongue cleaner
- Dental mirror

3. Instructions for Use

Product Description

1. Unit
2. Battery compartment cover
3. Stainless steel head
4. Silicone head
5. ON/OFF button
6. Tongue cleaner
7. Light
8. Dental mirror

BATTERY INSTALLATION

1. Turn the battery compartment cover (2) anti-clockwise and remove it.
 2. Insert the battery in accordance with the polarity marks on the unit and the cover (2).
 3. Replace the battery compartment cover (2) by aligning the part of the cover where there is an arrow with the part of the unit with a symbol of an open padlock. Press the cover (2) and turn it clockwise to the symbol of the closed padlock.
- Note: It is important for the battery compartment cover (2) to be closed properly so that contact is made and the product works. Ensure it is properly closed.

INSTRUCTIONS FOR USE

- The Starlyf® Ultrasonic is used to clean the front and back of teeth and between teeth.
- Choose the head that works best for you. First try the stainless steel head (3). If you don't like the feeling, try using the silicone head (4).
- Press the ON/OFF button (5) once to switch on the light (7).
- Press the ON/OFF button (5) again and the unit will vibrate while the light (7) stays on.
- If you press the button again, the unit will switch off.
- When the unit is vibrating, clean the teeth by simply pointing at and pressing on the area to be treated.
- Pull on the head to remove it and push down on the head to attach it.
- Use the mirror (8) to see all the areas to be treated.
- We recommend rinsing your mouth out after each cleaning.
- To clean your tongue, attach the tongue cleaner (6) to the unit.
- First clean the centre of your tongue using the green part of the cleaner with the outer crest.
- Scrape your tongue with the tongue cleaner from the centre to the tip.
- After each scrape, rinse the tooth cleaner with clean water and repeat scraping several times for maximum effect.
- Clean the sides of the tongue with the white side of the tongue cleaner. Repeat the operation for maximum effect.
- Clean your tongue before brushing your teeth.
- For reasons of hygiene, replace the head every three months.
- In the first days of use, the tongue cleaner may cause nausea in some cases. This sensation will go away with time and you'll be able to reach the back part of the tongue without difficulty.

4. MAINTENANCE AND CLEANING

IMPORTANT: Do not immerse the product in water.

- Ensure the unit is in the OFF position before cleaning it.
- First remove the heads from the unit (1) and then clean them with water or alcohol wipes.

QUALITY WARRANTY

This product is covered by a warranty against manufacturing defects subject to the time periods stipulated by the legislation in force in each country. This warranty does not cover damage resulting from inadequate use, negligent commercial use, abnormal wear and tear, accidents or improper handling.

Australia & New Zealand only

This product has the benefit of certain Consumer Guarantees. These are prescribed by the Australian Consumer Law & the New Zealand Consumer Guarantees Act 1993 both of which provide protection for consumers. There is no express warranty for this product in Australia or New Zealand. The above paragraph refers to other countries

Disposal of materials

A symbol of a crossed-out wheeled bin means you should find out about and follow local regulations about disposing of this kind of product.

Do not dispose of this product as you would other household waste. Dispose of this device in accordance with the corresponding local regulations. Electrical and electronic devices contain hazardous substances that can have harmful effects on the environment and/or human health and should be recycled properly.

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