FR - COUPE COMPRIMÉ 3EN1 - 127006 - VERSION 2025-10-07

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GB - PILL BOX WITH CUTTER 3-IN-1 - 127006

Pill box with four compartments for organizing and storing medication or vitamins. Includes pill cutter and mortar for easily halving/crushing medication.

Instructions for use:

The pill box has two different sides:

- one side with four compartments for storing medication.
- one side with a pill cutter and mortar.

Using the cutter/mortar:

Open the transparent lid of the cutter and place a tablet in the mold. Press down firmly on the transparent lid to split the tablet.

Open the transparent lid of the mortar and place a tablet on the grid-shaped surface. Press down firmly on the transparent lid to grind the tablet.

NOTE: Clean the pill cutter after each use (wash by hand). Only suitable for hard capsules and hard tablets.

- Handle with care as the blades are very sharp.
- The appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Keep out of reach of children.

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