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Folding Table RFM

USER MANUAL

Carefully read the whole of the user instructions before you use the folding table RFM for the first time. Keep the user instructions in a safe place so that they are available for later use.

Assembly/Instructions

1. Take all parts out of the cardboard packaging and check them for completeness and damage.
2. Included are: 1 x Table top, 1 x Bracket for adjusting and fixing the table top, 1 x u-shaped tube with push buttons at both ends (frame), 2 x l-shaped tube (frame), 1 x cross brace (frame), 2 x screw, one instruction manual

3. Insert an I-shaped tube to each end of the u-shaped tube and press the push buttons. Make sure that the snap fasteners engage through the openings at the same height.
4. Fix the cross brace between the two I-shaped tubes with the two screws. This serves to stabilize the frame. The frame should stand upright.

5. Hook the bracket, which serves to fix the table top, in the two openings on the left and on the right side of the frame.

6. First press the table top hanging on the standing frame.

7. Fold the table top up and hook the bracket underneath the table top for the desired angle. Use the push buttons to set the appropriate height.

Cleaning

You should regularly clean and disinfect the folding table RFM, which is intended for reuse. Use a ph-neutral soap solution and warm water. Dry all elements with a dry cloth. The disinfection can be carried out using a commercially available agent in accordance with the manufacturer's instructions. Do not use any abrasive or corrosive cleaning agents or high pressure or steam cleaners!

Technical Data/Dimensions and Weight

Weight: total 3,1kg, table top: 1 kg

Height: 53 - 74 cm

Dimensions table top: 52 x 40 cm

Width: 51 cm

Depth: 45 cm

Pipe diameter: 19 mm, 21 mm

Maximum load of the table top: 5 kg

Table top tilt angle: 0°, 15°, 30

Diameter of the drink holder: above 8 cm, below 7,5 cm

Material: steel, polypropylene

Guarantee

The guarantee period is 24 months and relates to all product defects which can be verifiably shown to be attributable to manufacturing or material defects.

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Rehaforum Medical GmbH, Daimlerstr. 12a, 25337 Elmshorn