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# Activon Tube™ GBR

*Activon honey in a 25g tube*

***Product description***
Pure sterile Manuka honey in a tube, effective against a broad spectrum of bacteria. Will maintain a moist wound healing environment and effectively de-slough and debride wounds. Will effectively reduce the odours associated with malodorous wounds.

***Indications***
May be used on all wound types including superficial wounds, partial or full thickness wounds, sloughy wounds, necrotic wounds e.g. pressure ulcers, leg ulcers, surgical wounds, burns, graft sites, fungating tumours.

***Directions for use***
May be applied directly to a wound or used with other Activon dressings.To apply, twist off the cap and squeeze tube, ensuring coverage of the wound area, replace cap and dispose of tube. To increase the loading of other Activon dressings apply to the dressing before applying the dressing to wound. Cover with a secondary dressing of choice.

***Frequency of change***
Depending on wound exudate levels, any surrounding interstitial fluid, oedema and dressing regimen Activon may remain for up to 7 days.Will need to be replenished as it is absorbed or diluted by wound exudate.To maintain an effective level re-apply as necessary. Due to its natural colour Activon is usually apparent within a wound; colour change may indicate absence of Activon.To remove excess Activon e.g. for wound inspection, irrigate wound with saline solution.

***Contra-Indications***
Product contains glucose. Monitor the levels of patients with diabetes. Do not use if allergic to bee venom. Pain may be experienced when applied, depending on sensitivity of the wound it may be necessary to consider an appropriate level of analgesic. If pain continues unabated, discontinue use and irrigate the wound with saline solution.

***Disposal***
Dispose as clinical waste.

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